

Ch(7): pH-metric titration of weak acids with strong bases and of weak bases with strong acids.

Objectives

- Write the equation of the reaction of weak acid with strong base.
 - Plot the curve $\text{pH} = f(V_b)$.
 - Determine the equivalence point from the titration curve.
 - Calculate the concentration C_a .
- Describe the curve $\text{pH} = f(V_b)$ and indicate the remarkable points.

I-pH-metric titration of weak acids (HA) with strong bases (OH⁻)

Consider a solution of weak acid (HA) having unknown concentration $C_a = ?$

A volume $V_{a \text{ taken}} = 10 \text{ mL}$ from this solution is placed in a beaker and titrated with a solution of strong base (NaOH) of known concentration $C_b = 10^{-2} \text{ mol/L}$

The acid/base pair related to HA is : (HA/A⁻) having a given pK_a


The pair related to strong base (NaOH) is : (H₂O/OH⁻) : $pK_a = 14$

The equation of titration reaction is :





This reaction is spontaneous (gamma rule), complete and should be quantitative ($K_R > 10^4$) and unique.

The pH of the solution is measured after each mL of NaOH added from the buret by using the pH-meter and the titration curve obtained has the following shape :



Curve $\text{pH} = f(V_b)$



Calculation of concentration C_a

The volume $V_{b \text{ equivalence}}$ is determined from the curve $\text{pH} = f(V_b)$ by the two parallel tangents method. At equivalence, acc. to SR :

$$n(\text{HA})_{\text{taken}} = n(\text{NaOH})_{\text{added at equivalence}}$$

$$n(\text{HA})_{\text{taken}} = n(\text{OH}^-)_{\text{added at equivalence}}$$

$$C_a \times V_{a \text{ taken}} = C_b \times V_{b \text{ eq}}$$

$$C_a = \frac{C_b \times V_{b \text{ eq}}}{V_{a \text{ taken}}}$$

Justification of the value of pH at equivalence point E_2

At equivalence, the two reactants (HA) and (OH^-) are reacted completely and transformed into products.

The species present in the titration beaker at equivalence are : A^- , H_2O and Na^+ obtained from NaOH.

Na^+ : is spectator ion that has no effect on the pH.

H_2O : is neutral.

A^- : is a weak base that reacts with water in a reversible reaction to produce a **little amount of OH^-** ions that gives a **pH > 7** according to the following reaction:



Justification of the value of pH at the half- equivalence point E_1

At the half equivalence : $V_{b \text{ added}} = \frac{V_{beq}}{2}$

At this point, the half of concentration of the acid HA is reacted and the other half is transformed into A^- , thus $[HA]_{\text{remained}} = [A^-]_{\text{formed}}$

Since the two species are present in the same solution, therefore the pH of the solution is determined by the Handerson relation :

$$pH = pK_a (HA/A^-) + \log\left(\frac{[A^-]}{[HA]}\right) = pK_a + \log(1) = pK_a$$

$$pH_{1/2 \text{ eq}} = pK_a (HA/A^-)$$

Description of the curve $pH = f(V_{b1})$

This curve is an ascending curve that consists of **4 parts** and **two inflection points** :

Part AB : The pH increases relatively rapidly.

Part BC : The pH increases slightly forming a plateau. This part consist of inflection point E_1 which represents the half equivalence.

Part CD : The pH increases rapidly and the curve shows a jump of pH. This part consists of the inflection point E_2 which represents the equivalence point.

Part DE : The pH increases slightly again and the curve tends to reach a limit of pH.

II-Titration of weak base (B^-) with strong acid (H_3O^+)

The solution of **weak base (B^-)** of unknown concentration C_b is placed in the beaker and titrated by a **strong acid (HCl or HNO_3)** of known concentration C_a .

The volume of the base taken is $V_{b \text{ taken}} = 10 \text{ mL}$ is placed in the beaker and the pH-meter measures the pH after each mL of the acid added from the buret.

The acid/ base pair related to B^- is: **(BH/B^-)** having a given pK_a .

The pair related to the strong acid (HCl or HNO_3) is: **(H_3O^+/H_2O)** having a $pK_a = 0$


The equation of titration reaction is :





$$K_R = 10^{(pK_a - 0)} = 10^{pK_a} > 10^4$$

This reaction is **spontaneous** (gamma rule), **complete** and should be **quantitative** ($K_R > 10^4$) and **unique**.

The titration curve obtained has the following shape :



Curve : $\text{pH} = f(V_a)$



Calculation of concentration C_b

The volume $V_{a \text{ equivalence}}$ is determined from the curve $\text{pH} = f(V_a)$ by using the two parallel tangents method.

$$n(\text{B}^-)_{\text{taken}} = n(\text{HCl})_{\text{added at equivalence}}$$

$$n(\text{B}^-)_{\text{taken}} = n(\text{H}_3\text{O}^+)_{\text{added}}$$

$$C_b \times V_{b \text{ taken}} = C_a \times V_{a \text{ equivalence}}$$

$$C_b = \frac{C_a \times V_{a \text{ eq}}}{V_{b \text{ taken}}}$$

Justification of the value of pH at the equivalence point E_2

At equivalence, the two reactants (B^-) and (H_3O^+) are reacted completely and transformed into products.

The species present in the titration beaker at equivalence are : **BH**, **H_2O** and **Cl^-** obtained from HCl.

Cl^- : is a spectator ion that has no effect on the pH.

H_2O : is neutral.

BH : is a weak acid that reacts with water in a reversible reaction to produce a **little amount of H_3O^+** ions that gives a **pH < 7** according to the following reaction:



Justification of the value of pH at the half equivalence point E₁

At the half equivalence : $V_{a \text{ added}} = \frac{V_{a \text{ eq}}}{2}$

At this point, the half of concentration of the base B⁻ is reacted and the other half is transformed into BH ,thus $[B^-]_{\text{remained}} = [BH]_{\text{formed}}$

Since the two species are present in the same solution, therefore the pH of the solution is determined by the Handerson relation :

$$\text{pH} = \text{pK}_a (\text{BH}/\text{B}^-) + \log \left(\frac{[B^-]}{[BH]} \right) = \text{pK}_a + \log (1) = \text{pK}_a$$

$$\text{pH}_{1/2 \text{ eq}} = \text{pK}_a (\text{BH}/\text{B}^-)$$

Description of the curve pH = f(V_a)

This curve is an descending curve that consists of **4 parts** and **two inflection points** :

Part AB : The pH decreases relatively rapidly.

Part BC : The pH decreases slightly forming a plateau. This part consist of inflection point E₁ which represents the half equivalence.

Part CD : The pH decreases rapidly and the curve shows a jump of pH. This part consists of the inflection point E₂ which represents the equivalence point.

Part DE : The pH decreases slightly again and the curve tends to reach a limit of pH.

III- Acid/base indicators

An acid base indicator is an **acid/ base pair (HA/A^-)** having a **given pK_a** where the conjugate **acid HA** has a color different than that of its conjugate **base A^-** .

For example : **HA** has a **yellow color** and **A^-** has a **blue color**.

Three cases are possible :

First case

If the indicator is placed in a solution having a **$\text{pH} < \text{pK}_a - 1$** , thus **$\text{HA}$** **predominates** in the solution and the color of the solution **becomes yellow**.

Second case

If the indicator is placed in a solution having a **$\text{pH} > \text{pK}_a + 1$** , thus **$\text{A}^-$** **predominates** in the solution and the color of the solution **becomes blue**.

Third case

If the indicator is placed in a solution having a pH that ranges between $\text{pK}_a - 1$ and $\text{pK}_a + 1$ (**$\text{pK}_a - 1 < \text{pH} < \text{pK}_a + 1$**), thus **$[\text{HA}] \approx [\text{A}^-]$** and the color of solution is a mixture of yellow and blue (**green**).

Examples of acid/base indicators

Indicators	pKa	Color of the acid species (HA)	pH change range	Color of the base species (A ⁻)
Methyl orange	3.6	Red	3.1 – 4.4 Orange	Yellow
Methyl red	5.1	Red	4.2 – 6.2 Orange	Yellow
Bromothymol blue	6.8	Yellow	6 – 7.6 Green	Blue
Phenolphthalein	9	Colorless	8.2 – 10 Pink	Purple